
aseptomag

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Table of contents

I.	Safety Instructions.....	4
II.	Valve Indication.....	5
II.I	General Information.....	5
II.II	Indications	6
III.	Installation	7
II.III	Welding Components with Pipe Connection	7
II.IV	Welding Components with Flange Connection	8
II.V	Welding Instructions "Flange Connection"	9
II.VI	Welding Tools.....	11
	Welding Tool type "O"	11
	Welding Tool type "S"	11
	Welding Tool to double chamber bottom seat valve S.....	12
	Welding Tool to double chamber bottom seat valve O.....	12
	Welding Tool to shut-off bottom-seat valve S	12
	Welding Tool to double-disc bottom seat valve S	12
II.VII	Welding Finalization.....	13
IV.	Compressed Air Supply	14
V.	Cleaning	14
VI.	Sterilization	14
VII.	Maintenance.....	15
	Planned Maintenance	15
	Precautionary Maintenance	16
	Maintenance After Leakage.....	17
	Lubrication	17
VIII.	Recommended Use of Bellows-Valves	18
IX.	Technical Data	18
	Operating Data.....	18
	Materials	19
	Surfaces	19
	Materials Certifications	20
X.	Trouble-shooting.....	20

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Always pay attention to the yellow safety data sheet!

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Valves and process systems represent a considerable hazard and can cause injuries. Work with extreme care when valves are connected to the air supply and/or process actions are taking place.

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Valves and process systems which are equipped with electrical equipment represent a considerable hazard and can cause injuries. Work with extreme care when installations and valves are connected to the power supply.

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Pay attention to the technical data! Only trained and qualified personnel should operate or maintain the valves.

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Never touch the piping, connections, or valves when they are hot! Never open the valve until the installation is empty and in a depressurized condition. Observe Lock-out/Tag-out.

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Each Aseptomag component is identified by the following reference numbers.

The numbers correlate to each other. For this reason, each number can be used to identify a component.

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Item Number (EDP-Nr.)

Each part and each assembly is identified by an article number with nine numerals. This is the primacy reference for a clear identification of a component.

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Secondary Reference

The Secondary Reference is composed of five groups. The first two letters stand for the valve type, the second group defines the nominal diameter and the three remaining groups describe the drawing numbers of the main components.

- ✓ AV = Aseptic Shut-off Valve
- ✓ 50 = Nominal diameter
- ✓ V50-1001 = Valve Housing
- ✓ V50-2200 = Internal Assembly
- ✓ V50-3000 = Pneumatic Actuator

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Drawing Number

The drawing number is a combination of two groups. The numerals in front of the hyphen refer to the nominal diameter and the remaining four numbers describe the part in further detail. The following defines the drawing number organisation for the main component groups of a valve:

- ✓ Vxx-0xxx = Entire Valve
- ✓ Vxx-1xxx = Valve Housing
- ✓ Vxx-2xxx = Internal Assemblies
- ✓ Vxx-3xxx = Actuators
- ✓ Vxx-4xxx = Feedback Units

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
An Aseptomag Valve, regardless of the valve type, generally contains four main components.

- ✓ Valve Housing
- ✓ Internal Assembly
- ✓ Actuator
- ✓ Clamp

Each of these main components are clearly marked with reference numbers so to be easily identified.

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Valve Housing

Drawing number Materials data for the housing identification are engraved in the same place. The serial number is located in a separate area (e.g. 054810).

Clamp

The clamp can be identified by checking the drawing number clearly cast-in the stainless steel.

Internal Assembly

Drawing and serial number of the Internal Assembly are engraved on the Valve Cover.

Pneumatic Actuator

Brief description, drawing, article and serial numbers of the Actuator are engraved on the cylinder.

Verifying Sticker

The serial number of the entire valve is printed on the round Verifying Sticker.

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11.111 Welding Components with Pipe Connection

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The pneumatic actuator, internal assembly and body seals must be removed before welding the valve into an installation. Please refer to the manual Mounting and Dismounting of Aseptomag Valves.

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Valve Activation

Warning: Before loosening the clamp, valves with direction of activation NC (normally closed) must be placed into the "open" position (Danger of injury from released spring force.)

Machining of Pipe Ends

Saw pipe ends flat and square, deburr and degrease before welding.

Weld in

Component housings must be welded in de-energized and free from distortion.

Welding Method

Apply welding method TIG (tungsten inert gas) or GTAW (gas tungsten arc welding) with butt weld I-joint according DIN8532. Hand or orbital weld.

Backing gas

Apply backing gas

Tacking

Tack the pieces under backing and shielding gas on 4 positions allocated around the periphery (Assure the backing and shielding gas supplies).

Welding Ends

Adjoined welding ends must not be separated by gap, because the exhausted backing gas reduces the corrosion resistance of the welded joint and the pipe.

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11.IV Welding Components with Flange Connection

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Tank Drill

Draw up the tank drill according to the valve specification with a tolerance from $d^{+0.5\text{mm}}$. With the accordant plate thickness a V-joint must be prepared.

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Weld in Equipment

For heat rejection reasons Aseptomag advises to use weld in equipment in general. They can be borrowed for free.

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Tacking

After mounting the suitable welding tool, tack the flange over the cross. This assures a distortion-free welding.

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
Welding Method

TIG (tungsten inert gas) or GTAW (gas tungsten arc welding). Tank flanges must be welded by pulse welding.

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11.V Welding Instructions "Flange Connection"

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
Tack Welding Process

Tack weld as illustrated in diagram above in a star pattern (1-2; 3-4; 5-6; 7-8). This will assure a distortion-free weld.

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Inner / Outer Welding Process

Seam type: V-Seam with bottom run

- ✓ Bottom run without ad material, top run with ad material (amount for top run according to material strength).
- ✓ Weld must be executed with a welding device with forming gas connected.
- ✓ Weld must be executed with the "step-back" operation.
- ✓ After every weld, cool seam with air or water
- ✓ Depending on welding tool type, the device must be turned before welding the reverse side.

IMPORTANT:

- ✓ Must be executed only with pulse welding
- ✓ Minimize current to reduce heat

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

Curing Interior

As long as the forming went well, a weld curing is not necessary in the interior of pipe connections.

Curing External Area

Curing method for external area: bating, brushing, grinding, polishing

Curing Tank Interior

According to the technical requirements, the welds at the tank interior will be ground or polished.

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After mounting it must be assured that the pipes are free of material residues such as saw dust, needless sealing materials, etc.

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To assure the optimal functionality of the components, Aseptomag recommends that the valve be installed so that the stream flows against the valve disk.

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For the valve to function correctly it needs to be operated with 6 bar (87psi) filtered (min. 0.5 µm) oil free compressed air.

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The valve is designed for CIP cleaning (Cleaning In Place).

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Recommended Cleaning Solution

- Chlorine-free water
- Nitric acid 0.8-1.2 %
- Caustic soda 1.2-2 %

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Recommended Cleaning Speed in the Valve
Minimum 2 m/sec

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The valve is designed for SIP sterilization (Sterilization In Place)

Sterilisation mit:

- Heisswasser max. 150 °C
- Dampf max. 150 °C
- Chemikalien (H₂O₂, Oxonia etc.)

Sterilization with:

- Hot water max. 150 °C
- Steam max. 150 °C
- Chemicals (H₂O₂, Oxonia, etc.)

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ASEPTOMAG components do not need any maintenance under regular working conditions. The only exceptions are wear parts such as seals and internal assemblies.
Always keep a sufficient stock of replacement wearing parts (internal assemblies and seals) on-hand.

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Always handle valve parts with the greatest care.

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Before mounting or dismantling always read the operating and maintenance instructions.

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After maintenance work always check the valve for correct operation by connecting the actuator to the compressed air supply and by opening and closing the valve several times.

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- ✓ Periodic check every 6 months
- ✓ Keep a maintenance log for the valves
- ✓ Use statistical analysis to plan maintenance intervals

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- ✓ Depending on the medium, check seal once a week.

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Gehäuse/Housing

- ✓ Check the interior of the housing for any damage.

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- ✓ Before mounting, grease the seals with PARA-LIQ GTE 703 (O/Art.No. 9999.50013)

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VIII. Recommended Use of Bellows-Valves

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Aseptomag bellows valves have been successfully placed in aseptic and hygienic installations since 1968. The regular integration and design is based on a maximum (mathematical) static operating pressure of 12 bar (174psi.) Aseptomag is able to design bellows valves up to 80 bar (1160psi.) Please contact our technicians regarding such a specialized design and integration.

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It is important to consult the technical data during mounting, operation, and maintenance. The following data references standard valve applications.

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Operating Data

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6 bar (andere gem. Spezifikation)

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Absolutes Vakuum

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Mindestens 6 bar (NW125/150 mindestens 8 bar)

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NC = federschliessend/luftöffnend (normally closed)

NO = luftschliessend/federöffnend (normally open)

L/L = luftöffnend/luftschliessend

Maximum Operating Pressure ΔP

6 bar (others acc. specification)

Minimum Operating Pressure

Absolute vacuum

Control Air Pressure Pneumatic Actuator

Minimum 6 bar (DN125/150 minimum 8 bar)

Direction of Action

NC = spring closing/air opening (normally closed)

NO = air closing/spring opening (normally open)

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Valve Seat Seal Materials (in contact with product)

- TEPASEP® (max. operating temp. 160 °C)
- PTFE
- Silicone
- EPDM
- Viton™

Housing Seal Materials (in contact with product)

- Silicone (max. operating temperature 160 °C)
- G-blau
- EPDM
- FEP-Silicone
- Viton™
- HNBR

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Further sealing material options upon request.

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Inside Area in contact with product

Surface roughness $R_a \leq 0.8 \mu\text{m}$ (standard)

Outer Surface

Metal blank fine turned and/or polished

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If desired, inside areas in contact with product can be electropolished/passivated or ground. These methods will bring up surface roughness to about $R_a \leq 0.6 / 0.4 \mu\text{m}$.

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All of the problem examples used below are known today. For problems not listed, please contact our customer service.

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